



Delicot®

THE UNEQUIVOCAL LEADER IN PLIABILITY




*"I have been using **Delicot**® sponges for almost twenty years now. There is no doubt in my practice that the soft consistency and pliability of these sponges has fundamentally shaped the way that I do surgery. Protection of neural and vascular structures just cannot be done as effectively or safely when using thicker, stiffer sponges. Therefore, I have come to rely on **Delicot**® sponges as a critical element to intracranial surgery and spinal cord microsurgery."*


*J.D. Day, M.D.
Professor and Chairman, Department of Neurosurgery*

Delicot® in a Right Side Vestibular Schwannoma Resection Retrosigmoid Approach (R)

Delicot®



BEING USED AS A SPACER TO LET OUT CSF



PROTECTING CRANIAL NERVE VII


63x00
2 x 20mm
Delicot®

SEE THE VIDEO:



Delicot® is the only product focused on delicate microsurgery procedures. With an ultrathin profile, 0.33mm thick, **Delicot**® is the most pliable patty on the market today. This allows for precise positioning, molding, and pliability for instrument support and tissue protection. **Delicot**® is capable of absorbing 8.7 times its dry weight in fluid while still demonstrating the best filterability via suction. All **Delicot**® patties are x-ray detectable and equipped with a locator string.

Delicot® Sizes



CONTACT US FOR
A FREE SAMPLE:



FOR MORE INFORMATION VISIT WWW.AMERICANSURGICAL.COM/DELICOT.