

Telfa®*

American Surgical Company's Precut Telfa® Patties Save Operating Room Time and Reduce the Potential for Retained Foreign Objects.

Telfa® dressings are commonly used by physicians looking for absorbency and protection of the underlying tissue while reducing the potential for adhesion, especially during prolonged use. Standard Telfa® sheets can address some of these needs, but their limited size offerings can lead to the practice of intra-operative cutting of larger Telfa® sheets in order to meet surgical requirements. The lack of x-ray markers and locator strings provides few options for the clinician to reduce the potential risk of the cut dressings being left behind after the procedure.



Until now...

The precision cut Telfa® neurosurgical patties and strips from ASC address these issues and provide the following advantages:

X-ray Detectable - *All ASC precut Telfa® patties contain an X-ray detectable strip for increased traceability and safety during extended procedures.*

Locator Strings - *Each ASC precut Telfa® pattie is equipped with a locator string and comes packed ten patties per pouch, incorporating a simple and efficient counting card.*

Available in 12 Sizes - *ASC's precut Telfa® patties and strips offer an extensive range of popular sizing options eliminating the need to spend time and resources manually cutting strips in the OR.*

Trusted Performance - *Manufactured to the same high standards as all ASC neurological patties, from a material you know and trust for low-adhesion.*

For more information or to request a free sample, please contact our Customer Service at 800-343-0060, Option 2, or visit our website at www.AmericanSurgical.com.