



2025 PRODUCT CATALOG





NEURO PATTIES

RAY-COT®	2
DELICOT®	4
TAPERED DELICOT®	5
UNIQCOT™	6
TELFA®	7
AMERICOT®	8

SINGLE-USE DEVICES

BIPOLAR FORCEPS	10
CRANIAL PERFORATORS	12


OTHER SPECIALTIES

COTTON BALLS	13
OPHTHALMIC PLEDGETS	13

RAY-COT® Rayon Patties

THE SAME TRUSTED PROPERTIES OF TRADITIONAL RAYON PATTIES, BUT STRONGER AND MORE PLIABLE

- Single x-ray detectable barium line allows for superior pliability and softness
- Sewn in locator string for increased safety
- Packaged in quantities of 200 per box (20 pouches of 10)



#60-02

6 mm x 6 mm
¼ in x ¼ in



#60-04

6 mm x 76 mm
¼ in x 3 in



#60-05

10 mm x 76 mm
¾ in x 3 in



#60-06

13 mm x 13 mm
½ in x ½ in



#60-07

13 mm x 25 mm
½ in x 1 in



#60-08

13 mm x 38 mm
½ in x 1½ in




#60-09

13 mm x 76 mm
½ in x 3 in



#60-10

19 mm x 19 mm
¾ in x ¾ in



#60-12

25 mm x 25 mm
1 in x 1 in




#60-13

25 mm x 76 mm
1 in x 3 in



#60-14

38 mm x 38 mm
1½ in x 1½ in



#60-15

76 mm x 76 mm
3 in x 3 in



#60-16

6 mm x 38 mm
¼ in x 1½ in



#60-22

38 mm x 76 mm
1½ in x 3 in





#60-29

13 mm x 51 mm
½ in x 2 in

RAY-COT® Rayon Strips

THE SAME TRUSTED PROPERTIES OF TRADITIONAL RAYON STRIPS,
BUT STRONGER AND MORE PLIABLE


- Single x-ray detectable barium line allows for superior pliability and softness
- Strips come without a locator string
- Packaged in quantities of 200 per box (20 pouches of 10)



DELICOT[®] Rayon Patties

THE ONLY PRODUCT FOCUSED ON MICROSURGICAL PROCEDURES AND THE LEADER IN PLIABILITY

- Ultra-thin 0.33 mm profile
- Single x-ray detectable barium line allows for superior pliability and softness
- Sewn in locator string for increased safety
- Packaged in quantities of 200 per box (20 pouches of 10)*



#63x00

2 mm x 20 mm
 $\frac{5}{16}$ in x $2\frac{3}{32}$ in



#63-01

5 mm x 8 mm
 $\frac{3}{16}$ in x $\frac{5}{16}$ in



#63-02

5 mm x 13 mm
 $\frac{3}{16}$ in x $\frac{1}{2}$ in



#63-03

10 mm x 13 mm
 $\frac{3}{8}$ in x $\frac{1}{2}$ in



#63-04

15 mm x 20 mm
 $\frac{9}{16}$ in x $1\frac{3}{16}$ in



#63-05

4 mm x 50 mm
 $\frac{5}{32}$ in x $1\frac{15}{16}$ in



#63-06

6 mm x 50 mm
 $\frac{1}{4}$ in x $1\frac{15}{16}$ in




#63-07

9 mm x 50 mm
 $\frac{11}{32}$ in x $1\frac{15}{16}$ in



#63-08

13 mm x 60 mm
 $\frac{1}{2}$ in x $2\frac{3}{8}$ in



#63-09

20 mm x 60 mm
 $1\frac{3}{16}$ in x $2\frac{3}{8}$ in



#63-10


30 mm x 70 mm
 $1\frac{3}{16}$ in x $2\frac{3}{4}$ in

*#63x00 Packaged in quantities of 100 per box (10 pouches of 10)

DELICOT® Tapered Rayon Patties


TAPERED DELICOT® MAXIMIZES PROTECTION IN THE MOST CONFINED OPERATIVE CORRIDORS.

- Tapered shape for increased performance
- Single x-ray detectable barium line allows for superior pliability and softness
- Sewn in locator string for increased safety
- Packaged in quantities of 200 per box (20 pouches of 10)




#63-07T

4 mm x 9 mm x 48 mm




#63-09T

8 mm x 19 mm x 64 mm



#63-12T

8 mm x 19 mm x 76 mm




#63-08T

6 mm x 19 mm x 51 mm




#63-13T

6 mm x 13 mm x 9 mm



#63-14T

19 mm x 25 mm x 25 mm




#63-15T

19 mm x 25 mm x 76 mm

THE SAME TRUSTED PROPERTIES OF TRADITIONAL RAYON PATTIES, BUT THINNER FOR BETTER VISIBILITY

- Uniqcot™ is the golden mean between Ray-Cot® and Delicot®
- Single x-ray detectable barium line allows for superior pliability and softness
- Sewn in locator string for increased safety
- Packaged in quantities of 200 per box (20 pouches of 10)



#67-02

6 mm x 6 mm
¼ in x ¼ in



#67-04

6 mm x 76 mm
¼ in x 3 in



#67-06

13 mm x 13 mm
½ in x ½ in



#67-07

13 mm x 25 mm
½ in x 1 in



#67-08

13 mm x 38 mm
½ in x 1½ in



#67-09

13 mm x 76 mm
½ in x 3 in



#67-10

19 mm x 19 mm
¾ in x ¾ in



#67-12

25 mm x 25 mm
1 in x 1 in



#67-13

25 mm x 76 mm
1 in x 3 in



#67-29

13 mm x 51 mm
½ in x 2 in

UNIQCOT™ Strips
No locator string



#67-81

6 mm x 152 mm
¼ in x 6 in



#67-82

13 mm x 152 mm
½ in x 6 in



#67-84

25 mm x 152 mm
1 in x 6 in

TELFA® Surgical Patties & Strips

A NATURAL/POLYMER HYBRID NEUROSURGICAL PATTIE FOR BOTH REDUCED ADHESION AND HIGH ABSORBENCY

- Telfa® outer layer surrounding an absorbent cotton layer
- Remarkably low tissue adhesion with high absorbency
- Single x-ray detectable barium line allows for superior pliability and softness
- Sewn in locator string for increased safety
- Packaged in quantities of 200 per box (20 pouches of 10)



Telfa® is a registered trademark of KPR U.S., LLC



#80-04
6 mm x 76 mm
¼ in x 3 in



#80-05
10 mm x 76 mm
⅜ in x 3 in



#80-06
13 mm x 13 mm
½ in x ½ in




#80-08
13 mm x 38 mm
½ in x 1½ in



#80-09
13 mm x 76 mm
½ in x 3 in



#80-11
19 mm x 76 mm
¾ in x 3 in



#80-12
25 mm x 25 mm
1 in x 1 in



#80-13
25 mm x 76 mm
1 in x 3 in



#80-14
76 mm x 76 mm
3 in x 3 in

TELFA® Strips
No locator string



#81-06
13 mm x 152 mm
½ in x 6 in




#81-13
25 mm x 152 mm
1 in x 6 in

AMERICOT® Cotton Patties

AMERICOT® IS CAPABLE OF HOLDING GREATER THAN 11 TIMES ITS DRY WEIGHT IN FLUID WHEN FULLY SATURATED

- 100% American Cotton
- Single x-ray detectable barium line allows for superior pliability and softness
- Sewn in locator string for increased safety
- Packaged in quantities of 200 per box (20 pouches of 10)



#20-06s
13 mm x 38 mm
½ in x 1½ in



#20-05s
13 mm x 25 mm
½ in x 1 in



#20-04s
13 mm x 13 mm
½ in x ½ in



#20-01s
6 mm x 6 mm
¼ in x ¼ in



#20-19s
3 mm x 3 mm
⅛ in x ⅛ in



#20-07s
13 mm x 76 mm
½ in x 3 in



#20-22s
10 mm x 10 mm
¾ in x ¾ in



#20-08s
19 mm x 19 mm
¾ in x ¾ in



#20-09s
25 mm x 25 mm
1 in x 1 in



#20-10s
25 mm x 76 mm
1 in x 3 in



#20-13s
19 mm x 76 mm
¾ in x 3 in



#20-12s
76 mm x 76 mm
3 in x 3 in



#20-14s
10 mm x 19 mm
¾ in x ¾ in



#20-27s
19 mm x 51 mm
¾ in x 2 in



#20-25s
10 mm x 25 mm
¾ in x 1 in



#20-26s
10 mm x 51 mm
¾ in x 2 in



#20-32s
13 mm x 51 mm
½ in x 2 in




#20-30s
6 mm x 38 mm
¼ in x 1½ in

AMERICOT[®] Cotton Strips

AMERICOT[®] IS CAPABLE OF HOLDING GREATER THAN 11 TIMES ITS DRY WEIGHT IN FLUID WHEN FULLY SATURATED

- 100% American Cotton
- Single x-ray detectable barium line allows for superior pliability and softness
- Strips come without a locator string
- Packaged in quantities of 200 per box (20 pouches of 10)



#21-01

6 mm x 152 mm

¼ in x 6 in



#21-02

13 mm x 152 mm

½ in x 6 in



#21-04

25 mm x 152 mm


1 in x 6 in

ATLAS™ CHOICE

Ultra Non-stick Bipolar Forceps

THE ATLAS™ CHOICE DISPOSABLE NON-STICK BIPOLAR FORCEPS FEATURE TIPS CREATED SPECIFICALLY FOR MICROSURGERY (MICRO TIPS) AND TUMOR AND SPINE (MACRO TIPS) PROCEDURES.

- Solid silver precision tips
- Stainless steel body for improved durability
- Low-profile slim tines
- Dual irrigation options
- Packaged in quantities of 1 per box



Length

Tip Size and Product Codes

8"

0.5 mm Micro Tip
Non-Irrigating: BIP-ATMI-04

1.0 mm Micro Tip
Non-Irrigating: BIP-ATMI-05
Irrigating: BIP-ATMI-05-IR

1.5 mm Macro Tip
Non-Irrigating: BIP-ATMA-03
Irrigating: BIP-ATMA-03-IR


2.0 mm Macro Tip
Non-Irrigating: BIP-ATMA-04
Irrigating: BIP-ATMA-04-IR

9"

0.5 mm Micro Tip
Non-Irrigating: BIP-ATMI-07
Irrigating: BIP-ATMI-07-IR


1.0 mm Micro Tip
Non-Irrigating: BIP-ATMI-08
Irrigating: BIP-ATMI-08-IR

1.5 mm Macro Tip
Non-Irrigating: BIP-ATMI-09
Irrigating: BIP-ATMI-09-IR



Superior Tip Functionality

Our solid silver tips create a superior non-stick effect. Perfect tip alignment provides balance and control for seamless dissection without a scissoring effect.



Advanced Heat Transport

Advanced heat transport technology efficiently regulates tine and tip temperature to minimize any tissue damage or charring.

MERIDIAN Cranial Perforators

MERIDIAN SINGLE-USE CRANIAL PERFORATOR FOR BURR HOLE TREPANATION FROM ADEOR MEDICAL AG.

- Unmatched sharpness
- Universal Hudson fitting
- Durable cutting efficiency
- Start & stop when needed
- Packaged in quantities of 1 per box



MERIDIAN CLASSIC
PER C6-9-S
6.0 mm / 9.0 mm (inner / outer)

- > For use on skulls as thin as 2.0 mm
- > Alternative to Integra® / Codman® 26-1223



MERIDIAN CLASSIC
PER C8-11-S
8.0 mm / 11.0 mm (inner / outer)

- > For use on skulls as thin as 2.0 mm
- > Alternative to Integra® / Codman® 26-1222



MERIDIAN CLASSIC
PER C11-14-S
11.0 mm / 14.0 mm (inner / outer)

- > For use on skulls as thin as 2.0 mm
- > Alternative to Integra® / Codman® 26-1221



MERIDIAN ALPHA
PER 7-11-S
7.0 mm / 11.0 mm (inner / outer)

- > For use on skulls as thin as 3.0 mm
- > Alternative to AcraCut® 200-243



MERIDIAN ALPHA
PER 11-14-S
11.0 mm / 14.0 mm (inner / outer)

- > For use on skulls as thin as 3.0 mm
- > Alternative to AcraCut® 200-271 or 200-241



MERIDIAN ALPHA
PER 7-11-SP
7.0 mm / 11.0 mm (inner / outer)

- > For use on skulls as thin as 1.0 mm
- > Alternative to AcraCut® 200-253




MERIDIAN ALPHA
PER 11-14-SP
11.0 mm / 14.0 mm (inner / outer)

- > For use on skulls as thin as 1.0 mm
- > Alternative to AcraCut® 200-251

*Only available in the United States

X-RAY DETECTABLE Cotton Balls


- 100% American Cotton
- Each premium strung cotton ball contains an x-ray detectable marker for increased safety and accountability
- Double strung locator string allow for easier traceability post-procedure
- Packaged in quantities of 100 per box (20 pouches of 5)



X-Large

#30-05

32 mm
1 1/4 in
(Strung)



Large

#30-04

25 mm
1 in
(Strung)



Medium

#30-03

19 mm
3/4 in
(Strung)



Small

#30-00

6 mm
1/4 in
(Strung)



X-Small

#30-07

3 mm
1/8 in
(Strung)



X-Small

#31-07

3 mm
1/8 in
(Non-Strung)



Small

#31-00

6 mm
1/4 in
(Non-Strung)

COMFORTPACK™ Ophthalmic Pledgets

THE COMFORTPACK™ EYE PLEDGETS FOR OPHTHALMOLOGY

- Absorbs anesthetic solutions for precise and effective numbing
- Eliminates waste and increases cost savings
- Sterile double strung cotton pledgets are a multi-use product for the clinic and operating room

#110-01

ComfortPack™ Kit

- > Two 5/32" (4 mm) double strung pledgets
- > One dispensing cup and swab
- > Packaged in quantities of 60 (30 pouches of 2)

#111-01

Ophthalmic
Pledgets 5/32"

- > Five 5/32" (4 mm) double strung pledgets only per pouch
- > Packaged in quantities of 100 (20 pouches of 5)

#111-02

Ophthalmic
Pledgets 7/32"

- > Five 7/32" (5.5 mm) double strung pledgets only per pouch
- > Packaged in quantities of 100 (20 pouches of 5)





Scan the QR Code to
Request a Sample Online

Order Now

1-800-343-0060 • 1-781-592-7200

CustomerService@AmericanSurgical.com

www.AmericanSurgical.com

